

# Lab (2)

## Bacterial Plasmid Extraction Method in the Laboratory

This laboratory outlines the procedure for extracting bacterial plasmids in a laboratory setting. Plasmids are small, circular DNA molecules found in bacteria that can replicate independently of chromosomal DNA. They are widely used in molecular biology for cloning, gene expression, and genetic engineering. The following steps provide a detailed methodology for plasmid extraction, ensuring high yield and purity.

### Materials Required

- Bacterial culture containing plasmid
- LB broth or appropriate growth medium
- Centrifuge
- Microcentrifuge tubes
- Lysis buffer (e.g., alkaline lysis buffer)
- Neutralization buffer
- Wash buffer (e.g., ethanol or isopropanol)
- Elution buffer (e.g., Tris-EDTA or water)
- Pipettes and tips
- Ice
- Vortex mixer
- Spectrophotometer (optional for quantification)

### Procedure

#### 1. Bacterial Culture Preparation

1. **Inoculation:** Start by inoculating a single colony of the bacterial strain containing the plasmid into 5-10 mL of LB broth supplemented with the appropriate antibiotic.
2. **Incubation:** Incubate the culture overnight at 37°C with shaking (200-250 rpm) to ensure optimal growth.

## 2. Cell Harvesting

1. **Centrifugation:** Transfer the overnight culture into a sterile microcentrifuge tube and centrifuge at 4,000-6,000 x g for 10 minutes at 4°C. This will pellet the bacterial cells.
2. **Supernatant Removal:** Carefully discard the supernatant without disturbing the pellet.

## 3. Cell Lysis

1. **Resuspension:** Resuspend the bacterial pellet in 200-300 µL of lysis buffer. Ensure the cells are fully resuspended by vortexing gently.
2. **Lysis:** Add 200 µL of alkaline lysis buffer to the resuspended cells and mix gently by inverting the tube several times. Allow the mixture to incubate at room temperature for 5 minutes. This step will lyse the cells and release the plasmid DNA.

## 4. Neutralization

1. **Neutralization:** Add 150-200 µL of neutralization buffer to the lysate. Mix gently by inverting the tube. This step will neutralize the alkaline conditions and precipitate cellular debris.
2. **Centrifugation:** Centrifuge the mixture at 12,000-14,000 x g for 10 minutes at room temperature. The plasmid DNA will remain in the supernatant, while cellular debris will form a pellet.

## 5. Plasmid Precipitation

1. **Supernatant Transfer:** Carefully transfer the supernatant containing the plasmid DNA to a new microcentrifuge tube.
2. **Precipitation:** Add an equal volume of ethanol or isopropanol to the supernatant and mix gently. Incubate the mixture at -20°C for at least 30 minutes to precipitate the plasmid DNA.
3. **Centrifugation:** Centrifuge at 12,000-14,000 x g for 10 minutes at 4°C. The plasmid DNA will form a pellet.

## 6. Washing the Plasmid DNA

1. **Wash:** Discard the supernatant and wash the pellet with 70% ethanol. Centrifuge again at 12,000-14,000 x g for 5 minutes.
2. **Drying:** Carefully remove the ethanol and allow the pellet to air dry for 5-10 minutes. Avoid over-drying, as this can make the DNA difficult to resuspend.

## 7. Resuspension of Plasmid DNA

1. **Resuspension:** Resuspend the dried plasmid DNA pellet in 30-50  $\mu\text{L}$  of Tris-EDTA or sterile water. Gently pipette up and down to ensure complete dissolution.
2. **Storage:** Store the plasmid DNA at  $-20^{\circ}\text{C}$  for long-term storage or at  $4^{\circ}\text{C}$  for short-term use.

## Quality Control

To confirm the successful extraction of plasmid DNA, it is advisable to perform a quality check using agarose gel electrophoresis. Additionally, measuring the concentration and purity of the plasmid DNA using a spectrophotometer can provide further validation.

The extraction of bacterial plasmids is a fundamental technique in molecular biology that allows researchers to manipulate genetic material for various applications. Following the outlined procedure will yield high-quality plasmid DNA suitable for cloning, transformation, and other downstream applications. Proper handling and storage of plasmid DNA are crucial for maintaining its integrity and functionality in subsequent experiments.